


CryoPro[®] Probes & Tips





Sharp pointed contact probe. Gold plated.
E.g. small angiomas, molluscum contagiosum.

1 mm contact probe. Gold plated.
E.g. small warts, angiomas, lesions near the eye, trichiasis.



2 mm contact probe. Gold plated.
E.g. small flat warts, lesions near the eye, trichiasis.

3 mm contact probe. Gold plated.
E.g. flat warts, lesions near the eye, lentiginos.



4 mm contact probe. Gold plated.
E.g. flat warts, lesions near the eye, lentiginos.

5 mm contact probe. Gold plated.
E.g. flat warts, lesions near the eye, lentiginos.




6 mm contact probe. Gold plated.
E.g. flat warts, lesions near the eye, lentiginos.

8 mm contact probe. Gold plated.
E.g. flat warts, lesions near the eye, lentiginos.



10 mm contact probe. Gold plated
E.g. squamous cell/ basal cell carcinoma, hemorrhoids.

15 mm contact probe. Gold plated
E.g. squamous cell/ basal cell carcinoma, hemorrhoids.



20 mm contact probe. Gold plated.
E.g. squamous cell/ basal cell carcinoma.

30 mm contact probe. Gold plated.
E.g. squamous cell/ basal cell carcinoma.



Bent spray.
E.g. for lesions difficult to reach, perianal, rectal or vaginal condyloma.

Straight spray.
E.g. for lesions difficult to reach, perianal, rectal or vaginal condyloma.



Acne soft spray tip.
For superficial desquamation, highly vaporized spray.

Luer lock adapter.
Allows attachment of accessories with luer lock fitting - e.g. disposable needles.



CORTEX TECHNOLOGY

Smedevaenget 10 9560 Hadsund Denmark

Tel.: +45 9857 4100 Fax: +45 9857 2223 cortex@cortex.dk www.cortex.dk
German office: Tel.: +49 2385 922090 Fax: +49 2385 922091 cortexde@aol.com